

Improved Protein Arrays for Quantitative Systems Analysis of the Dynamics of Signaling Pathway Interactions

Chin-Rang Yang

Staff Scientist, National Heart, Lung and Blood Institute, National Institute of Health

Summary

Astronauts and workers in nuclear plants who repeatedly exposed to low doses of ionizing radiation (IR, <10 cGy) are likely to incur specific changes in signal transduction and gene expression in various tissues of their body. Remarkable advances in high throughput genomics and proteomics technologies enable researchers to broaden their focus from examining single gene/protein kinetics to better understanding global gene/protein expression profiling and biological pathway analyses, namely Systems Biology. An ultimate goal of systems biology is to develop dynamic mathematical models of interacting biological systems capable of simulating living systems in a computer. This Glue Grant is to complement Dr. Boothman's existing DOE grant (No. DE-FG02-06ER64186) entitled "The IGF1/IGF-1R-MAPK-Secretory Clusterin (sCLU) Pathway: Mediator of a Low Dose IR-Inducible Bystander Effect" to develop sensitive and quantitative proteomic technology that suitable for low dose radiobiology researches.

An improved version of quantitative protein array platform utilizing linear Quantum dot signaling for systematically measuring protein levels and phosphorylation states for systems biology modeling is presented. The signals are amplified by a confocal laser Quantum dot scanner resulting in ~1000-fold more sensitivity than traditional Western blots and show the good linearity that is impossible for the signals of HRP-amplification. Therefore this improved protein array technology is suitable to detect weak responses of low dose radiation. Software is developed to facilitate the quantitative readout of signaling network activities. Kinetics of EGFRvIII mutant signaling was analyzed to quantify cross-talks between EGFR and other signaling pathways.

Keywords: protein array, quantum dot, quantitative analysis, signaling pathway, cross-talk, image analysis, glioma, systems biology

Background

The emerging field of Systems Biology focuses on the most daunting challenges in biology and medicine. It is an attempt to understand how all parts of the cell - genes, proteins, and other molecules - work in concert to create complex living organisms and analyzing how entire biological systems function, both in health and in sickness. System biology always relates big amount of data, which makes high throughput technologies are crucial. The rapid advance of high throughput technologies has enabled scientists to broaden their research from detailed investigation of a few selected genes/proteins to global gene/protein expression profiles and network analysis. Among the network analysis, cellular signal transduction networks play an important role in regulating cellular processes, such as proliferation, cell growth and death. Proteins are the work-horses that carry out these functions. Therefore, it is crucial to capture the dynamics of protein kinases and post-translational regulations within cellular signal transduction networks for understanding how the signaling pathways are operated in healthy versus disease conditions.

Reverse phase protein lysate array (RPPA), originally introduced by Drs L. Liotta and E. Petricoin [1], is designed for measuring protein expression in a large number of biological samples quantitatively. RPPAs are essentially high density arrays of micro protein lysate dot-blot spots spotted in series of dilutions that provide quantitative protein expression measurements on many samples in a single array. Arrays are probed with a primary antibody followed by a species-specific secondary antibody similar to the Western blot. The detection signal comes from the tag on the secondary antibody. A range of detection tags have been developed including colorimetric, fluorescent, near-infrared, and Quantum dot (Qdot) assays [2-6]. RPPA has been applied to protein monitoring for biomarker discovery and/or signal transduction proteins in response to various biological stimuli or chemical treatments [7-10]. However, to use RPPA as a quantification assay is a real challenge, because the linear signals, the foundation of quantification, are difficult to be obtained by using most detection methods.

When combined with the advanced quantum dot (Qdots), a nano-metal fluorophore with bright and linear signals and no photo-bleaching effect observed in organic fluorophores, RPPA becomes a sensitive quantification assay. Here, we present an enhanced version of the protein array platform for sensitive, reproducible transductional and quantitative cellular signal transduction network measurements. The cell lysis buffer is optimized for RPPA printing and dissolves whole cell proteins without using urea. The thin-coated-nitrocellulose slide is chosen for strong protein binding and low fluorescence background. Most importantly, we eliminated the common enzyme (horseradish peroxidase, HRP) signal amplification step in other protocols [2-5]. We demonstrated the enzyme amplification step is not linear, resulting in nonlinearity signals. In contrast, a confocal laser Qdot scanner is utilized to amplify and maintain the signal linearity. To further reduce background fluorescence from nitrocellulose and increase signal/noise ratios, the confocal laser is focused right above the nitrocellulose coating. Integrated software is used to automatically analyze array images, qualify and quantify spots in series, and generate serial dilution curves to determine the relative protein levels and phosphorylation states in the samples.

To demonstrate the capacity of our platform to capture the dynamics of signaling responses, and determine the sensitivity to detect minute changes, glioma cancer cells expressing constitutively activated EGFRvIII mutant under tetracycline control were analyzed by protein arrays. The EGFRvIII mutant is a common oncogenic mutant co-expressed with wild-type EGFR in glioblastoma (GBM) [11]. EGFRvIII is unable to bind ligand and signals constitutively. Kinetics of signaling after conditional induction of EGFRvIII expression was analyzed to quantify the response. The dynamics of pathway interaction (i.e. cross-talk) between EGFR pathways and other signaling pathways were then captured.

Results and discussion

Understanding complex cellular systems will require the identification and analysis of each of its components and allow determination of how they function together and are regulated. A critical step in this process is to determine the biochemical activities of the proteins and how these activities themselves are controlled and modified by other proteins. Traditionally, the biochemical activities of proteins have been elucidated by studying single molecules, one experiment at a time. This process is not optimal, as it is slow and labor intensive. To obtain a global view of molecular events instead of individual molecules, we utilize the protein array approach to monitor the molecular network response.

The advantages of the enhanced Qdot-RPPA platform are presented here.

Linear Dynamic [12] for Quantization To test the linearity of Qdot signal, we produced cell lysates with artificial gradient of p53 protein by mixing null p53 cell lysate (human lung cancer H358 cells with null p53 gene) and p53 mutant cell lysate (human lung cancer H2009 cells with overexpressing mutant p53) in different proportions. Figure 1A is the array image amplified by confocal Laser. Figure 1B is the array image amplified by HRP enzyme. As shown in Figure 1C, R^2 value of the Laser-amplified linear regression curve is close to 1 indicating Qdot-RPPA can distinguish at least 25% change of protein levels within the linear range, verse the signals of HRP-amplification showing the sigmoid response indicating lost of the linearity.

Sensitivity In Figure 1D, using purified Akt protein, Qdot-RPPA, like HRP-RPPA[, can detect as low as 0.1 pg compared to 0.1 ng of the detection limit using traditional Western blot in normal conditions, at least 1000-fold more sensitive. Moreover, the signal linearity over serial dilutions makes Qdot-RPPA a reliable tool for quantification (Figure 2).

Specificity Qdot-RPPA can detect specific kinase activities with validated phospho-specific antibodies. Qdot-RPPA detects phospho-AKT (pAKT ser473) specifically in the kinase- and phosphatase- inhibitor treatments (Figure 2A). Total AKT is not altered. Thus, with a pair of validated total and phospho-specific antibodies, Qdot-RPPA can be used to monitor functional status of a given kinase under different treatments or diseases.

Reproducibility The signals of Qdot are reliable over time. The same samples were spotted and hybridized with Qdot at the different time in duplicate. R^2 values are about 0.95.

Linearity Tests for Phospho-specific Antibodies and Caspase In Figure 2, a series of commercial available negative and positive controls for indicated phospho-specific and caspase 3 antibodies were mixed proportionally to generate artificial gradients of target proteins: 2A. pAKT; 2B. pERK; 2C. pGSK3; 2D. pNFkB; 2E. pp38; 2F. cleaved caspase 3. R^2 values of these Laser-amplified linear regression curves are near 1 assuring that Qdot-RPPA can distinguish at least 25% change of protein functional levels while using validated , specific antibodies.

Array Image Processing and Statistic Readouts

Compared to commercially available DNA microarrays, protein arrays present additional challenges in image analysis. The variety of array formats, spot shapes, and intensity profiles makes it challenging to extract spot signals correctly. In addition, the different array substrates, printing mechanisms and protocols, staining/blocking processes, and broad applications result in various kinds of complex images. It is extremely difficult to develop one algorithm fit all applications. MicroVigeneTM is implemented based on the object-oriented technology and enable a robust software system integrating with multiple algorithms that is flexible, configurable as well as extensible to provide customized solutions, support any future needs, and adapt along with this emerging field through proper plug-ins [13].

In Figure 3, customized MicroVigeneTM microarray image analysis software has been developed for high throughput, automatic array image processing and quantitative readouts from serial dilutions of samples. The software implements the actual boundary algorithm for spot identification/segmentation that is resistant to the spot shift and image shift; the regional background algorithm for local non-uniform background correction and sensitive spot quantification.

Also, instead of generating multiple linear regression curves for data quantification over each series of serial dilutions, MicroVigeneTM implements the SuperCurve algorithm that using all spots within one array to form a sigmoid antigen-antibody binding kinetic curve (i.e. SuperCurve) (Figure 3A). The advantage of SuperCurve is the resistance to the experimental outliers or missing spots on the array compared to regular linear regression curves formed by only a few spots over a serial dilution. The SuperCurve is a consensus curve supported by all spots on the array. This process has been implemented into the automatic data analysis workflow after the spot identification to automatically generate a reaction curve for each antibody as well as the quantitative readouts from each dilution series.

The automatic dust finding and removal algorithm is also implemented to increase the accuracy of curve fitting and the intensity readouts. In Figure 3B, a small but intensified dust spot was automatically detected in the lowest dilution spot, and then visualized by 3D image manually. Without the dust removal, the spot was flagged as an outlier of the SuperCurve verse. However, after removing the dust signal, the spot fitted in the SuperCurve well (Figure 3C).

In summary, these unique features of MicroVigeneTM make it capable of handling shifted and noisy protein array images and enabling the hand-free batch process required for high throughput protein array image processing.

The Dynamics of EGFRvIII Signaling Network

Effective targeted treatment of glioblastoma (GBM) may require a coordinated inhibition of multiple signals. Much remains to be learned about functional interactions among signal transduction networks in GBM. For example, GBMs express multiple receptor tyrosine kinase (RTK) families and ligands. Increasing evidence suggest extensive cross-talking between RTK signaling networks that have functional implications for multi-target treatment. EGFRvIII is a common oncogenic mutant that is co-expressed with the wild-type EGFR in GBM. We applied Qdot-RPPA to investigate key kinases that are regulated by conditional inducible EGFRvIII in U87MG glioma cell line. The dynamics of pathway interaction (i.e. cross-talk) among canonical EGFR pathway, Akt, Src and JNK pathways after turning on the EGFR vIII were captured (Figure 4). Qdot-RPPA detected the total EGFR, pEGFR and pERK levels changes within one to six hours after adding tetracycline; Activation of Akt [14, 15], Src [16, 17] and JNK [18-20] happened at later time (20-24 hrs), and, interestingly, pSrc and pJNK were activated at much higher levels than the downstream pERK in the canonical pathway. The results suggest the selection of potential candidates for the future multi-target treatment in GBM, e.g. co-targeting EGFRvIII and Src/JNK kinases for cancer treatment [16, 17].

Conclusions

In conclusion, we have set up an enhanced quantitative proteomic technology-Qdot-RPPA platform by using RPPA, in conjunction with the emerging Quantum dot (Qdot) nanotechnology, confocal laser scanning, computer automated array image analysis and quantitative outputs. This platform have completely get rid of HRP-amplification, while have still kept the advantages of RPPA such as: minute sample per spot (~1nl); high sensitivity (~0.1 pg); high specificity with validated antibodies; high reproducibility; and large capacity of spots on one slide allowing investigating a large number of experimental conditions in parallel within one array.

The primary contribution of Qdot-RPPA platform is to easily obtain linear signals without photo-bleaching effect and large dynamic range for quantification from protein array. This technology offers us the capacity to quantitatively monitor the time series and dose responses of cellular response over signaling network after treatments or different disease stages, facilitates the complex functional analysis among different signaling pathways. We also found that thin coated (10 um) nitrocellulose slide with low fluorescence background yields high signal/noise ratio.

We have examined Qdot-RPPA platform by conducting the Dynamics of EGFRvIII Signaling Network. The results showed that the platform worked well. The ultimate goal of this platform is to apply high-throughput functional protein array data in building a mathematical model for simulating the biological network response *in silico*

and provide direction for more in-depth experimentation prediction and hypothesis generation.

Methods

Cell Lysate Preparation

p53 control cell lysates were prepared from human H358 (p53 null) lung cancer cells and H2009 (p53 mutant) cells. The rest of control cell extracts were purchased from Cell Signaling Technology: Akt (cat#9273), p38 MAP kinase (cat#9213), NF- κ B (cat#9243), and caspase-3 (cat#9663). Purified Akt protein was purchased from Invitrogen (cat#P2999). U87MG glioma cells were stably transfected with an EGFRvIII mutant using a tetracycline inducible system as described previously [21]. Cells were treated with tetracycline (1 μ g/ml), and then, harvested at 0, 6, 20, 24 and 48 hours, respectively. The following cell lysis protocol is optimized to harvest total cell proteins including hard-to-dissolved membrane proteins without using Urea. Cells were washed twice in cold phosphate buffered saline (PBS) and subsequently lysed in 250-350 μ L hot lysis buffer (2% SDS; 0.06 M Tris-Cl, pH 6.8; 5% Glycerol). Proteinase inhibitor (Sigma, cat#p8340) and phosphatase inhibitor cocktails (Santa Cruz Biotechnology, sc-45044 and sc-45045) and 2% β -mercaptoethanol were freshly added before use. The cellular lysates were boiled for 5 min on heat block, followed by 1 min vortexing, and then centrifuged at 13,000 rpm for 7 min at 4°C. Cells should be completely dissolved without much precipitates left at this step. Supernatants were transferred to new tubes and stored at -20°C. Protein concentration was measured using Bradford assay (Bio-Rad, cat#500-0006). All the samples were adjusted to 0.5 μ g/ μ L as the highest concentration on the arrays. However, to correct for protein loading, Sypro Ruby™ protein stain signals were used as described below.

Qdot-Reverse-Phase Protein Microarray (Qdot-RPPA)

Protein lysates were filtered through 96-well filter plate with 25 μ m pore membrane (Phenix Research Products, cat#MPF-009) by centrifuge to remove sticky aggregates, and serially diluted (1:2) 4 times using 1 x lysis buffer. Lysates were arrayed on ONCYTE® AVID nitrocellulose film slides (Grace Bio-Labs, cat#305170) using a SpotArray™24 Microarray printing system (PerkinElmer) with 55-60% humidity. Spots were separated with 350 μ m space in between. Approximately 1 nl of lysate per spot was arrayed using 4 spotting pins (TeleChem/ArrayIt, cat#946MP3). Slides were dried at RT for about 30 min, and then stored at 4 °C. Immuno-staining was performed within one week. We have tested various slides from different vendors. The AVID slides were chosen due to its high protein binding capacity with less protein lost, and thin-coated (10 micron) nitrocellulose film with less auto-fluorescence background, suitable for both Qdot™ and SyproRuby™ stains.

Laser-amplified Qdot-RPPA and Immuno-staining

The slides were placed into four-chamber plates (ISC Bioexpress, cat#T-2896-1) and incubated at room temperature in Re-blot Plus Mild Solution (CHEMICON, cat#2502) for no more than 7 minutes to relax protein structure. After the Re-blot was removed, the slides were washed 5 min three times in TBS-T buffer (2.42 g Tris-HCl, 16

g NaCl in 1 L dH₂O, 0.1% Tween-100, pH 7.6). The slides were incubated in Sea Block blocking buffer (Thermo, cat#37527) at 4°C overnight, and then blocked with Avidin and Biotin (Dako, Biotin blocking system cat#X0590) sequentially. Between steps, slides were washed with TBS-T buffer. After the blocking, the slides were incubated with primary antibodies (diluted in antibody diluent buffer, Dako, cat#S3022) at 4°C overnight. We used 61 validated primary antibodies. Next, slides were washed with TBS-T buffer, and incubated with biotinylated secondary antibodies (Vector, anti-rabbit BA-1000 or anti-mouse BA-9200 IgG) (1:5000) for 30min, then with Qdot 655–streptavidin conjugate (Invitrogen, cat#Q10121MP) for 30 min. Qdot were diluted 1:200 in boric acid buffer (50mM Borate, 2% BSA, pH 8.3). The slides were washed with TBS-T buffer 3 times, and MilliQ water once and then briefly spun at 2000 rpm for 5min to dry the slides. To detect Qdot signals, slides were scanned with full Laser power on a ProScanArray Microarray Scanner (PerkinElmer), and the Laser was focused 10 micron above the slide surface. The excitation wavelength was set at 488 nm and emission wavelength was set at 655 nm using a Qdot 655 filter. Images were saved in 16-bit TIFF format and the maximum signal intensity was 65535.

HRP-amplified Qdot-RPPA

Qdot-RPPA using the horseradish peroxidase (HRP) based system (Dako CSA system, cat#K1500) to amplify the Qdot signal was previously described [5, 6]. This protocol was performed only for the comparison shown in Figure 1B.

Automatic Array Image analysis and Quantification

Scanned TIFF images were batch-analyzed using MicroVigene™ software (VigeneTech inc. <http://www.vigenetech.com/>). In RPPA, cell lysates were printed on the arrays. Each array is probed with one specific primary antibody and scanned for one image. It is necessary to have computer-aid image analysis that can process array images automatically to increase the throughput of the number of antibodies one can apply. This work was done in collaboration with VigeneTech, Inc., a leader of developing automated image analysis technology. MicroVigene™ provided unique software customized to analyze the Qdot-RPPA platform. Unique features of the software provide accuracy, sensitivity, and reliable results of automation [13] including flexible grid and actual spot boundary algorithm to quantify spot signals accurately; dust removal algorithm to remove the contaminated signals; regional background algorithm for local non-uniform background correction and sensitive spot quantification. Ten thousand spots can be processed in less than one minute and, with the feature of hands-free batch processing, MicroVigene™ enables the high throughput protein array image analysis. 3D visualization of processed images is also available for manually quality assurance.

Each sample has five dilutions and was printed in triplicate on the array. Therefore, each dilution series has total of 15 data points to minimize errors and increase the confidence of curve-fitting. Instead of generating multiple linear regression curves for data quantification over each series of serial dilutions, MicroVigene implements the SuperCurve algorithm (a 4-parameter logistic-log model, i.e. parameters a-d shown in equation below) that uses all spots within one array to form a sigmoid antigen-antibody binding kinetic curve (i.e. SuperCurve) [22].

$$Y = a + \frac{(b-a)}{(1 + e^{(c \cdot (d - \ln(x)))})}$$

where x is the dilution factor and Y is the signal intensity. The signal readout of each dilution series is the intensity of EC50 from the fitted SuperCurve. The assumption is that the same antibody-antigen binding kinetics is taking place at each sample spot, even in the different samples, thus by taking all spots on an array to fit a common response curve can increase the confidence of the curve fitting.

Calibration for protein loading was based on total protein per spot. For estimation of total protein amounts, randomly selected arrays were stained with Sypro Ruby™ (Invitrogen, cat#S11791), and visualized on ProScanArray Scanner (PerkinElmer) with excitation wavelength at 450 nm and emission wavelength at 610 nm. The corrected values were calculated by dividing the EC50 readouts of the antibody to the corresponding EC50 readouts of the Sypro Ruby™ stain.

List of abbreviations

GBM: glioblastoma

HRP: horseradish peroxidase

Qdot: Quantum dot

RPPA: Reverse phase protein lysate array

RTK: receptor tyrosine kinase

References

1. Liotta LA, Espina V, Mehta AI, Calvert V, Rosenblatt K, Geho D, Munson PJ, Young L, Wulfkuhle J, Petricoin EF, 3rd: **Protein microarrays: meeting analytical challenges for clinical applications.** *Cancer Cell* 2003, 3:317-325.
2. Spurrier B, Ramalingam S, Nishizuka S: **Reverse-phase protein lysate microarrays for cell signaling analysis.** *Nat Protoc* 2008, 3:1796-1808.
3. Kornblau SM, Tibes R, Qiu YH, Chen W, Kantarjian HM, Andreeff M, Coombes KR, Mills GB: **Functional proteomic profiling of AML predicts response and survival.** *Blood* 2009, 113:154-164.
4. Brase JC, Mannsperger H, Frohlich H, Gade S, Schmidt C, Wiemann S, Beissbarth T, Schlomm T, Sultmann H, Korf U: **Increasing the sensitivity of reverse phase protein arrays by antibody-mediated signal amplification.** *Proteome Sci*, 8:36.
5. Shingyoji M, Gerion D, Pinkel D, Gray JW, Chen F: **Quantum dots-based reverse phase protein microarray.** *Talanta* 2005, 67:472-478.
6. Geho D, Lahar N, Gurnani P, Huebschman M, Herrmann P, Espina V, Shi A, Wulfkuhle J, Garner H, Petricoin E, 3rd, et al: **Pegylated, streptavidin-conjugated quantum dots are effective detection elements for reverse-phase protein microarrays.** *Bioconjug Chem* 2005, 16:559-566.
7. Sheehan KM, Calvert VS, Kay EW, Lu Y, Fishman D, Espina V, Aquino J, Speer R, Araujo R, Mills GB, et al: **Use of reverse phase protein microarrays and reference standard development for molecular network analysis of metastatic ovarian carcinoma.** *Mol Cell Proteomics* 2005, 4:346-355.
8. Sheehan KM, Gulmann C, Eichler GS, Weinstein JN, Barrett HL, Kay EW, Conroy RM, Liotta LA, Petricoin EF, 3rd: **Signal pathway profiling of epithelial and stromal compartments of colonic carcinoma reveals epithelial-mesenchymal transition.** *Oncogene* 2008, 27:323-331.
9. VanMeter AJ, Rodriguez AS, Bowman ED, Jen J, Harris CC, Deng J, Calvert VS, Silvestri A, Fredolini C, Chandhoke V, et al: **Laser capture microdissection and protein microarray analysis of human non-small cell lung cancer: differential epidermal growth factor receptor (EGFR) phosphorylation events associated with mutated EGFR compared with wild type.** *Mol Cell Proteomics* 2008, 7:1902-1924.
10. Wulfkuhle JD, Aquino JA, Calvert VS, Fishman DA, Coukos G, Liotta LA, Petricoin EF, 3rd: **Signal pathway profiling of ovarian cancer from human tissue specimens using reverse-phase protein microarrays.** *Proteomics* 2003, 3:2085-2090.
11. Hatanpaa KJ, Burma S, Zhao D, Habib AA: **Epidermal growth factor receptor in glioma: signal transduction, neuropathology, imaging, and radioresistance.** *Neoplasia* 2010, 12:675-684.
12. Heimberger AB, Suki D, Yang D, Shi W, Aldape K: **The natural history of EGFR and EGFRvIII in glioblastoma patients.** *J Transl Med* 2005, 3:38.
13. Ruan M: **Protein Microarray Image Analysis.** In *Functional Protein Microarrays in Drug Discovery*. 1 edition. Edited by Predki PF: CRC Press; 2006: 359

14. Li L, Dutra A, Pak E, Labrie JE, 3rd, Gerstein RM, Pandolfi PP, Recht LD, Ross AH: **EGFRvIII expression and PTEN loss synergistically induce chromosomal instability and glial tumors.** *Neuro Oncol* 2009, 11:9-21.
15. C VS, Babar SM, Song EJ, Oh E, Yoo YS: **Kinetic analysis of the MAPK and PI3K/Akt signaling pathways.** *Mol Cells* 2008, 25:397-406.
16. Andersen P, Villingshoj M, Poulsen HS, Stockhausen MT: **Improved response by co-targeting EGFR/EGFRvIII and Src family kinases in human cancer cells.** *Cancer Invest* 2009, 27:178-183.
17. Bromann PA, Korkaya H, Courtneidge SA: **The interplay between Src family kinases and receptor tyrosine kinases.** *Oncogene* 2004, 23:7957-7968.
18. Antonyak MA, Kenyon LC, Godwin AK, James DC, Emlet DR, Okamoto I, Tnani M, Holgado-Madruga M, Moscatello DK, Wong AJ: **Elevated JNK activation contributes to the pathogenesis of human brain tumors.** *Oncogene* 2002, 21:5038-5046.
19. Antonyak MA, Moscatello DK, Wong AJ: **Constitutive activation of c-Jun N-terminal kinase by a mutant epidermal growth factor receptor.** *J Biol Chem* 1998, 273:2817-2822.
20. Aikin R, Maysinger D, Rosenberg L: **Cross-talk between phosphatidylinositol 3-kinase/AKT and c-jun NH2-terminal kinase mediates survival of isolated human islets.** *Endocrinology* 2004, 145:4522-4531.
21. Ramnarain DB, Park S, Lee DY, Hatanpaa KJ, Scoggin SO, Otu H, Libermann TA, Raisanen JM, Ashfaq R, Wong ET, et al: **Differential gene expression analysis reveals generation of an autocrine loop by a mutant epidermal growth factor receptor in glioma cells.** *Cancer Res* 2006, 66:867-874.
22. Hu J, He X, Baggerly KA, Coombes KR, Hennessy BT, Mills GB: **Non-parametric quantification of protein lysate arrays.** *Bioinformatics* 2007, 23:1986-1994.

Figures




Figure 1. Linearity of Laser-amplified Qdot-RPPA for quantifying protein levels.

A. Human H358 (p53 null) lung cancer cell lysate was mixed with H2009 (p53 mutant) cell lysate proportionally (%) as indicated, then printed with seven serial dilutions in triplicate on protein array. The array was probed for total p53 levels and signals amplified by confocal laser Qdot scanner.

B. The same lysates in A., but signals amplified by HRP.

C. Linear regression curve of 1:8 dilution from Panel A is shown here. X axis is % of p53 mutant lysate in the mixed sample and Y axis is the corresponding Qdot intensity normalized with total protein loading. Black squares is the data from Panel A; red diamond is the data from Panel B. 25% change of p53 can be detected by Laser-amplified Qdot-RPPA as indicated by the black linear regression line. In contrast, linearity was lost with HRP-amplification shown as a red sigmoid curve.

D. Qdot-RPPA could detect as low as 0.1pg of purified AKT. From top to bottom are 10 serial dilutions. From left to right represent 10 replicates. The Signal /Noise ratios are from 30 to 150 as indicated.

Figure 2




Figure 2. The Linearity tests for Quantifying Cellular Signaling Kinases

A series of negative (-) and positive (+) controls for phosph-specific and caspase 3 antibodies were mixed proportionally to generate artificial gradient of target proteins for the linearity tests using laser-amplified Qdot-RPPA. Examples shown here are (A) pAKT, (B) pERK, (C) pGSK3, (D) pNFkB, (E) pp38, (F) cleaved caspase 3. pAKT+, pERK+ and pGSK+ are Jurkat cells treated with phosphatase inhibitor, calyculin A, to accumulate these kinases in phosphorylated states. pNFkB+ is HeLa cells stimulated with TNF (20ng/ml, 5min) to activate NFkB. p-p38+ is C6 glioma cells treated with anisomycin to activate p38. Cleaved caspase 3+ is cytochrome c treated Jurkat cell lysates to induce the cleavage of caspase 3. X axis is the gradient (25% intervals) of targeted protein. Y axis is the calibrated signal intensity after subtracting the background signal at the 0% gradient.

R^2 values of linear regression curves close to 1 indicate laser-amplified Qdot-RPPA can distinguish at least 25% change of protein levels within the linear range.

Figure 3




Figure 3. Automation of Qdot-RPPA Image and Data Analysis. Several features of MicroVigeneTM software are to facilitate the automatic batch process. **A.** Super Curve: Instead of generating multiple linear regression curves for each serial dilution, the Super Curve algorithm uses all spots on the same array to fit one antigen-antibody reaction curve (called Super Curve, the red line) for the quantitative readouts. The purple lines indicate error boundary. **B.** 2D and 3D views of the actual spot boundary algorithm for actual spot identification and dust finding (indicated by asterisk); the total intensity within the defined spot boundary (i.e. volume) was used for Super Curve fitting. **C.** Screen shots of the dust removal feature. Left: before removing the dust, the spot is flagged as an outlier (the red spot pointed by the arrow). Right: after dust removal, the spot fits the Super Curve.

Figure 4




Figure 4. Signaling network up-regulated by constitutively activated EGFRvIII mutation U87MG glioma cells were transfected with a conditional inducible tetracycline (tet-on) sytem. A EGFRvIII expressing line (row 2) as well as vector control (row 1) was treated with tetracycline (1 ug/ml) for 0, 1, 6, 20, 24 and 48 hours as indicated. Activities of signaling kinases were monitored by Qdot-RPPA. Each sample was serially diluted 1:1 to 1:16 for quantification and spotted in triplicate on the arrays, then probed with total and phospho-specific kinase antibodies. Relative fold changes to time 0 were quantified and illustrated with color under each corresponding array images. Black arrows indicate known canonical signaling pathways; Red dashed arrows indicate the cross-talk between Akt, Src and JNK pathways. p38 and NFκB stay inactive.