

ORIGINAL ARTICLE

Open Access



Point-of-care detection of *Neisseria gonorrhoeae* based on RPA-CRISPR/Cas12a

Qianrong Tu¹, Xiaoying Cao², Chao Ling¹, Lili Xiang³, Ping Yang^{1*} and Shifeng Huang^{1*} 

Abstract

Gonorrhoea, caused by *Neisseria gonorrhoeae* (*N. gonorrhoeae*), is a persistent global public health threat. The development of low-cost, point-of-care testing is crucial for gonorrhoea control, especially in regions with limited medical facilities. In this study, we integrated CRISPR/Cas12a reaction with recombinase polymerase amplification (RPA) to provide a simple and adaptable molecular detection method for *N. gonorrhoeae*. The RPA-Cas12a-based detection system developed in this study enables rapid detection of *N. gonorrhoeae* within 1 h without the use of specialized equipment. This method is highly specific for identifying *N. gonorrhoeae* without cross-reactivity with other prevalent pathogens. Furthermore, in the evaluation of 24 clinical samples, the detection system demonstrates a 100% concordance rate with traditional culture, which is being used clinically as a reference method. Overall, the RPA-Cas12a-based *N. gonorrhoeae* detection has the advantages of rapidity, portability, low-cost, no special equipment required, and strong operability, and has a high potential for application as a self-testing and point-of-care diagnosis, which is critical for the clinical management of gonorrhoea in developing countries lacking medical equipment.

Key points

- Rapid detection of *N. gonorrhoeae* based on CRISPR/Cas12a and RPA.
- Specific and sensitive detection.
- No special equipment promises applications for point-of-care diagnosis.

Keywords *Neisseria gonorrhoeae*, Recombinase polymerase amplification, CRISPR/Cas12a, Fluorescence, Lateral flow detection


*Correspondence:

Ping Yang
201903@hospital.cqmu.edu.cn
Shifeng Huang
sfhuang@hospital.cqmu.edu.cn

Full list of author information is available at the end of the article



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>.



Introduction

Gonorrhoea is a sexually transmitted infection (STI) caused by the bacterium *Neisseria gonorrhoeae* (*N. gonorrhoeae*), with purulent inflammation of the urogenital tract as the main clinical symptom (Lovett and Duncan 2018; Stevens and Criss 2018). Moreover, gonorrhoea is associated with infertility, adverse pregnancy and newborn outcomes (Stevens and Criss 2018; Vallely et al. 2021), and it facilitates HIV acquisition and transmission (Stenger et al. 2021). Infrequently, *N. gonorrhoeae* infection can spread systemically, causing septic arthritis, skin manifestations (Lovett and Duncan 2018), and septicemia (WHO 2016). Gonorrhoea is a global public health challenge. According to the World Health Organization (WHO), the annual incidence of gonorrhoea cases among adults (15–49 years old) was estimated at 82.4 million globally in 2020 (WHO 2021), with the highest prevalence in low-income countries (Rowley et al. 2019). With a rapidly rising incidence (Qin et al. 2019; Yin et al. 2018), gonorrhoea has become one of the most frequently reported infectious diseases in China (Yin et al. 2018). Nationwide, the detection rate of *N. gonorrhoeae* is quite low despite its considerable prevalence. Traditional methods for *N. gonorrhoeae* identification mainly include culture and direct microscopy. Culture is labor-intensive and time-consuming; direct microscopy necessitates extensive laboratory technician experience and is prone to misdiagnosis. The PCR technique necessitates sophisticated equipment that may not be readily available in hospitals with minimal resources. Therefore, it is crucial to create a sensitive, low-cost, and rapid point-of-care test for *N. gonorrhoeae*.

Recent advancements in clustered regularly interspaced short palindromic repeats (CRISPR)/associated proteins

(Cas) systems offer a new exciting and promising avenue for molecular diagnostics. On the basis of CRISPR/Cas systems, researchers have successfully developed a variety of sensitive, specific and rapid detection platforms for detecting various pathogens, such as SHERLOCK (Cas13a) (Gootenberg et al. 2017), DETECTR (Cas12a) (Chen et al. 2018), and Cas14-DETECTR (Harrington et al. 2018). Cas12a (Cpf1), an RNA-guided endonuclease, has both cis- and trans-cleavage DNase activities (Zetsche et al. 2015; Swarts 2019). After targeted recognition of DNA sequences complementary to CRISPR RNA (crRNA), a cis-cleavage activity is displayed by Cas12a, which cleaves double-stranded DNA (dsDNA) (Swarts 2019). Subsequently, the activated Cas12a displays collateral trans-cleavage activity, that is, nonspecific cleavage of nearby untargeted single-stranded DNA (ssDNA), an activity extensively explored for use in nucleic acid diagnostics (Swarts 2019). By introducing fluorophore and quencher-labeled or fluorophore and biotin-labeled ssDNA into the reaction system, CRISPR/Cas12a cleavage results can be detected by fluorescence readout or lateral flow test strips (Xiong et al. 2020).

The isothermal amplification feature of recombinase polymerase amplification (RPA) reduces the equipment needs and opens up new avenues for breaking past laboratory boundaries. In addition, RPA's tolerance to background DNA and certain (PCR) inhibitors such as haemoglobin, heparin, and urine (Li et al. 2018), which facilitates its clinical application, is of great significance.

The *porA* pseudogene was often selected as the target for *N. gonorrhoeae* detection due to its high specificity (Mangold et al. 2007). The *porA* gene of *N. meningitidis* encoded porin protein, while the *porA* pseudogene was not expressed in *N. gonorrhoeae*. As a result of the

absence of positive selection pressure, the pseudogene appears to be extremely stable and has remained nearly unchanged over time (Hjelmevoll et al. 2006). In addition, the *porA* gene/pseudogene was found only in two human-pathogenic *Neisseria* species, *N. meningitidis* and *N. gonorrhoeae* (Hjelmevoll et al. 2006). Importantly, the *porA* pseudogene of *N. gonorrhoeae* to be identified is sufficiently distinct from the *porA* gene of *N. meningitidis* (Hjelmevoll et al. 2006).

Based on the DETECTOR platform, this study has integrated CRISPR/Cas12a reaction with RPA to establish a new nucleic acid diagnostic method for *N. gonorrhoeae*. The RPA-Cas12a system enables rapid detection of *N. gonorrhoeae* without the need for specialized equipment, which is of great significance for the diagnosis and control of gonorrhea in developing countries lacking medical equipment.

Materials and methods

Reagents

Primers for RPA were synthesized by TsingKe Biological Technology (Beijing, China). Recombinant RNase inhibitor, crRNA, and RNase-free water were obtained from Takara Bio (Beijing, China). The ssDNA reporter was synthesized by Sangon Biotech (Shanghai, China). The TIANamp Bacteria DNA kit was purchased from TIANGEN (Beijing, China). The TwistAmp Basic Kit for RPA reaction was purchased from TwistDx (Cambridge, UK). Lba Cas12a (Cpf1) and NEBuffer r2.1 were purchased from New England Biolabs (Beijing, China). Milenia HybriDetect for lateral flow detection was purchased from Milenia Biotec (Gießen, Germany).

Bacterial strain and clinical samples

N. gonorrhoeae reference strain (ATCC 43,069) was donated by Professor Xiaobing Zhang. *N. gonorrhoeae* clinical isolates (n=8) were identified by Gram staining and mass spectrometry. 10 non-*N. gonorrhoeae* clinical isolates, including *Neisseria sicca* (*N. sicca*, n=2), *N. meningitidis* (n=2), *Ureaplasma urealyticum* (*U. urealyticum*, n=1), *Pseudomonas aeruginosa* (*P. aeruginosa*, n=1), *Enterobacter cloacae* (*E. cloacae*, n=1), *Klebsiella pneumoniae* (*K. pneumoniae*, n=1), *Enterococcus faecium* (*E. faecium*, n=1), and *Acinetobacter radioresistens* (*A. radioresistens*, n=1), were selected as negative controls. In addition, 24 clinical samples were collected from patients suspected of having gonorrhea, including cervical or vaginal discharge from women and urethral secretions from men. Sufficient clinical samples were manually separated into two equal halves, one for traditional culture (gold standard) and one for DNA extraction for molecular diagnosis using the RPA-Cas12a system. All the clinical isolates and samples described above were

obtained from the First Affiliated Hospital of Chongqing Medical University.

RPA primers, crRNA and ssDNA reporter design

All the available *porA* pseudogene sequences of *N. gonorrhoeae* (Accession Numbers AJ010732.1, AJ010733.1, AJ223449.1, AJ223448.1, AJ223447.1, and AJ223446.1) were downloaded from the NCBI database (www.ncbi.nlm.nih.gov) and aligned by Mega-X (www.megasoftware.net) to identify the conserved regions. To distinguish *N. gonorrhoeae* from *N. meningitidis*, we aligned the conserved sequences of the *porA* pseudogene of *N. gonorrhoeae* with the *porA* gene of 228 *N. meningitidis* strains (the NCBI accession numbers of *porA* gene of *N. meningitidis* were listed in Additional file 1: Table S1). RPA primers (Additional file 1: Table S2) were designed in accordance with the TwistDx instruction manual (www.twistdx.co.uk) to cover the conserved *porA* pseudogene region that differs significantly from the *porA* gene sequences of *N. meningitidis*. IDT OligoAnalyzer (www.idtdna.com) (Kersting et al. 2014) was used to examine the primer dimer and hairpin structures. The primers' specificities were verified by NCBI's Primer-BLAST (<https://www.ncbi.nlm.nih.gov/tools/primer-blast>). We designed crRNA (Additional file 1: Table S3) for Cas12a to recognize a 21-bp target sequence near the PAM site (TTTN, N=A/C/G) in the amplicon between RPA primer pairs, and checked the crRNA sequence using the BLAST algorithm. Fluorescent ssDNA reporter was labeled by fluorophore and quencher (Additional file 1: Table S4). Lateral flow ssDNA reporter was labeled by fluorophore and biotin (Additional file 1: Table S4).

Nucleic acid preparation

Bacterial genomic DNA was extracted with the TIANamp Bacteria DNA kit per the manufacturer's instructions. The DNA concentration ranged from 20 ng/μL to 100 ng/μL, and was stored at -20 °C before use after being determined using a NanoDrop 2000 (Thermo, United States). For sensitivity testing, a solution containing genomic DNA from a *N. gonorrhoeae* reference strain was diluted to various concentrations with RNase free water.

The RPA reaction

To achieve robust assay performance and a specific RPA reaction, a comparative evaluation of the specificity of multiple candidate RPA primer pairs in identifying *N. gonorrhoeae* was carried out. Briefly, six candidate primer pairs were used to amplify *N. gonorrhoeae* and non-*N. gonorrhoeae* strain sequences, and the optimal primer pairs were evaluated by nucleic acid electrophoresis or the CRISPR/Cas12a-fluorescent assay. Using

Table 1 The RPA reaction system

Component	Volume (μL)
Primer F (10 μM)	2.4
Primer R (10 μM)	2.4
Primer-free rehydration buffer	29.5
RNase free water	11.2
Template	2

Add the above reaction mixture to a tube containing lyophilized RPA enzyme. Then add 2.5 μL 280 mM MgOAc before incubation at 39 °C for 15–20 min

Table 2 The RPA-Cas12a assay system

Component	Volume (μL)
NEBuffer r2.1 (10 \times)	2
Cas12a (1 μM)	1
crRNA (1 μM)	1.2
Fluorescent ssDNA reporter (2 μM), or Lateral flow ssDNA reporter (5 μM)	1
RNase inhibitor	0.25
RNase free water	13.55
RPA products	1

a commercial RPA kit, the RPA reaction was carried out at 39 °C for 15–20 min (Table 1) according to the instructions.

RPA-Cas12a-fluorescent assay

Fluorophore and quencher-labeled ssDNA reporter was introduced into the Cas12a-fluorescent system. Following the RPA reaction, the RPA product was added as substrate into the Cas12a-fluorescent system. The RPA-Cas12a-fluorescent reaction was carried out at 37 °C in a final volume of 20 μL (as specified in Table 2). The fluorescence signal was examined by TECAN infinite 200Pro (TECAN, Switzerland) with an excitation wavelength of 492 nm, and an emission wavelength of 522 nm. The optimization of cleavage time within the range of 5 to 55 min was undertaken to attain desirable performance of the CRISPR/cas12a system.

RPA-Cas12a-LFD assay

Except for the ssDNA reporter, the RPA-Cas12a-lateral flow detection (LFD) system contained the same components as the RPA-Cas12a-fluorescent assay (Table 2). In addition, after CRISPR/Cas12a cleavage at 37 °C, the CRISPR/Cas12a reaction product was diluted proportionately in Hybri Detect Assay Buffer to meet the solution volume requirements for LFD. The strips were then inserted and incubated at room temperature for 2 min before being withdrawn and photographed. ImageJ (<https://imagej.nih.gov/ij/>) quantification and GraphPad

(www.graphpad.com) visualization were used to examine the band intensities of the test line of the strips. Several parameters of the RPA-Cas12a-LFD system were optimized, such as the ssDNA concentration (3.125 nmol/L, 6.25 nmol/L, 12.5 nmol/L, 25 nmol/L, 62.5 nmol/L, 125 nmol/L, 250 nmol/L), the volume ratio of Cas12a cleavage product to Hybri Detect Assay Buffer (1:1, 1:3). Commercial strips were used for lateral flow detection. The Milenia HybriDetect kit contains the Hybri Detect Assay Buffer and lateral flow strips.

Specificity and sensitivity of the RPA-Cas12a system

Using DNA templates from the *N. gonorrhoeae* reference strain (ATCC43069), eight *N. gonorrhoeae* clinical isolates and ten non-*N. gonorrhoeae* clinical isolates (including *N. sicca*, *N. meningitidis*, *U. urealyticum*, *P. aeruginosa*, *E. cloacae*, *K. pneumoniae*, *E. faecium*, *A. radioresistens*), the specificity of the RPA-Cas12a system was evaluated. Following a 20-min RPA reaction, the RPA products were added into both the Cas12a-fluorescent system and the Cas12a-LFD system.

The sensitivity of the RPA-Cas12a system was evaluated using serial dilutions of genomic DNA of the *N. gonorrhoeae* reference strain (500 pg/ μL , 50 pg/ μL , 5 pg/ μL , 500 fg/ μL , 50 fg/ μL , 5 fg/ μL). The RPA products were added into both the Cas12a-fluorescent system and the Cas12a-LFD system following a 20-min RPA reaction.

Clinical validation of the RPA-Cas12a system

The feasibility of the RPA-Cas12a system was evaluated using 24 clinical samples from individuals suspected of having gonorrhea, with traditional culture as a reference method. Sufficient clinical samples were manually separated into two equal halves, one for traditional culture, and the other for DNA extraction for molecular diagnosis using RPA-Cas12a-fluorescent and RPA-Cas12a-LFD assays. In other words, all clinical samples were evaluated using the aforementioned three methods, and the test results were compared with those of the culture method. Culturing was performed on Thayer-Martin (TM) medium in 5% CO₂ at 37 °C for 24–72 h (Su et al. 2011).

Results

Construction of specific RPA reaction

The specificities of the six RPA primer pairs used to pre-amplify the *porA* pseudogene of *N. gonorrhoeae* were initially validated using NCBI's Primer-BLAST. In the process of identifying the optimal primer pairs, we discovered an intriguing phenomenon. Candidate primer pairs were used to amplify *N. gonorrhoeae* and non-*N. gonorrhoeae* strain sequences, and the nucleic acid electrophoresis results indicated that the primer pair No. 3 and No. 6 were specific (Fig. 1a, b). However, when the




Fig. 1 Optimizations of the RPA primer pairs and the CRISPR/Cas12a cleavage time durations. **a** Evaluation of the specificity of the RPA primer pair No.3 by electrophoresis on a 2% agarose gel. 1, *N. gonorrhoeae* ATCC43069; 2, *E. faecium*; 3, *N. sicca*; 4, *N. sicca*; 5, *N. meningitidis*; 6, *N. meningitidis*; 7, *P. aeruginosa*; 8, *A. radioresistens*; 9, *E. cloacae*; 10, *K. pneumoniae*. **b** Evaluation of the specificity of the RPA primer pair No.6 by electrophoresis on a 2% agarose gel. 1–7, *N. gonorrhoeae* clinical isolates; 8, *A. radioresistens*; 9, *K. pneumoniae*; 10, *N. meningitidis*; 11, *N. meningitidis*. BP, base pair; M, marker. **c** Evaluation of the specificity of the RPA primer pairs by Cas12a detection system. RPA reactions were performed with 6 sets of primer pairs. **d** Optimizations of the CRISPR/Cas12a cleavage time durations. Ng, *N. gonorrhoeae* ATCC43069 (1 ng/ μ L); Nm, *N. meningitidis*. AU, arbitrary unit

RPA amplification product generated by primer pair No. 6 was added to the Cas12a detection system, the specificity of the detection system was compromised. Therefore, we abandoned the direct screening of specific primer pairs by electrophoresis and instead observed the specificity in conjunction with the Cas12a-fluorescent detection system to select the most suitable RPA primer pairs for the Cas12a detection system. As depicted in Fig. 1c, the primer pair No. 3 exhibited a high degree of specificity, which is advantageous for the Cas12a detection system.

Construction of the RPA-Cas12a-fluorescent assay

Under the condition of producing distinct fluorescence intensity, the shorter the CRISPR/Cas12a cleavage time,

the quicker the RPA-Cas12a-fluorescent assay. Therefore, the cleavage time was optimized between 5 and 55 min, and because the fluorescence was consistent at 30 min, 30 min was chosen as the optimal time. However, the fluorescence intensity of *N. gonorrhoeae* was readily distinguishable from that of the negative control strain *N. meningitidis* before reaching a plateau. Thus the duration of cleavage can be reduced to 15 to 30 min (Fig. 1d).

Construction of the RPA-Cas12a-LFD assay

Although the RPA-Cas12a-fluorescent assay was easy-to-implement, it was dependent on the fluorescent readout instrument, therefore, we developed the RPA-Cas12a-LFD assay to achieve the on-site diagnosis of *N. gonorrhoeae*, with the advantages of portability, low cost and

no need of specialized equipment. To reduce false-positive results in the RPA-Cas12a-LFD assay, we diluted ssDNA into different concentrations. ImageJ quantification and GraphPad visualization of the band intensities of the test line of lateral flow strips revealed that 125 nmol/L ssDNA produced the lowest band intensity for the test line of negative samples (Fig. 2a, b). In addition, the volume ratio between the Cas12a reaction mixture and the Hybri Detect Assay Buffer was also optimized. The Cas12a reaction mixture was diluted with the Hybri Detect Assay Buffer in different proportions (1:1, 1:3), and the findings demonstrated that a ratio of 1:1 can make the positive band the most prominent (Fig. 2c).

Specificity and LoD of RPA-Cas12a system

DNA extracts from *N. gonorrhoeae* reference strain (ATCC43069), *N. gonorrhoeae* clinical isolates (n=8) and non-*N. gonorrhoeae* clinical isolates (*N. sicca*, *N. meningitidis*, *U. urealyticum*, *P. aeruginosa*, *E. cloacae*, *K. pneumoniae*, *E. faecium*, *A. radioresistens*) were employed as templates to evaluate the specificity of the RPA-Cas12a system. Both the RPA-Cas12a-fluorescent and RPA-Cas12a-LFD assays can accurately identify *N. gonorrhoeae* without cross-reactivity with non-*N. gonorrhoeae* (Fig. 3a, b, d), indicating high specificity.

To determine the limit of detection (LoD) of the RPA-Cas12a system, serially diluted *N. gonorrhoeae* ATCC43069 DNA extracts were used as templates (500 pg/ μ L, 50 pg/ μ L, 5 pg/ μ L, 500 fg/ μ L, 50 fg/ μ L, 5 fg/ μ L). The results showed that the LoD of the RPA-Cas12a-fluorescent assay and RPA-Cas12a-LFD assay were both as low as 5 pg/ μ L (Fig. 3c, e), demonstrating high sensitivity.

Clinical validation of the RPA-Cas12a system

In addition, the feasibility of the RPA-Cas12a-based *N. gonorrhoeae* detection system was evaluated using clinical samples from 24 individuals suspected of having gonorrhea. When taking traditional culture as a reference method, both RPA-Cas12a-fluorescent assay and RPA-Cas12a-LFD assay demonstrated a 100 percent concordance rate (Fig. 4a, b; Table 3). Our method promises clinical applications for the identification of *N. gonorrhoeae* due to its high feasibility.

Discussion

Current clinical methods for detecting *N. gonorrhoeae*, whether the traditional culture, direct microscopy, or the emerging PCR method, all require specialized equipments and skilled technicians, limiting the convenience and adaptability of the detection. The RPA-Cas12a-based *N. gonorrhoeae* detection system developed in this study provides rapid detection of *N. gonorrhoeae* without the need for specialized equipment and has significant potential for application as a point-of-care test, particularly the RPA-Cas12a-LFD assay. Notably, *N. gonorrhoeae* and *N. meningitidis* share a high degree of homology (Vigué and Eyre-Walker 2019), making it difficult to distinguish between the two pathogens during detection. The RPA-Cas12a-based *N. gonorrhoeae* detection method targeting the *porA* pseudogene had no cross-reactivity with other pathogens, and accurately distinguished *N. gonorrhoeae* from *N. meningitidis*. Therefore, its value for clinical *N. gonorrhoeae* detection is significant. In addition to its great specificity, the RPA-Cas12a system demonstrated a high sensitivity with a limit of detection (LoD) of 5 pg/ μ L. Besides, the practicability of the detection system was evaluated using 24 clinical samples, showing a 100%




Fig. 2 Construction of the RPA-Cas12a-LFD assay **a, b** The optimization of the concentration of the lateral flow detection ssDNA reporter (3.125 nmol/L, 6.25 nmol/L, 12.5 nmol/L, 25 nmol/L, 62.5 nmol/L, 125 nmol/L, 250 nmol/L). C, control line; T, test line. Band intensity, intensity of the test line was measured by ImageJ. **c** The optimization of the ratio of Cas12a reaction mixture to Hybri Detect Assay Buffer. Cas12a reaction mixture: Hybri Detect Assay Buffer (1:1, 1:3)




Fig. 3 Specificity and LoD of RPA-Cas12a system **a, b** Specificity of RPA-Cas12a-fluorescent assay, with reaction lasting for 30 min. **c** The LoD of RPA-Cas12a-fluorescent assay, with reaction lasting for 30 min. **d** Specificity of RPA-Cas12a-LFD assay, with reaction lasting for 30 min. **e** The LoD of RPA-Cas12a-LFD assay, with reaction lasting for 30 min. Data points represent at least four biologically independent experiments. Error bars indicate the mean \pm SD. AU, arbitrary unit. The asterisk indicates a significant difference compare with the control of no targets by one-way ANOVA analysis (* $P < 0.05$, **** $P < 0.0001$). ns, not statistically significant. Ng, *N. gonorrhoeae* ATCC43069; Nm. 1, *N. meningitidis*. 1; Nm. 2, *N. meningitidis*. 2; Ns. 1, *N. sicca*. 1; Ns. 2, *N. sicca*. 2; Pae, *P. aeruginosa*

concordance with the traditional culture method. Our method promises clinical applications for the identification of *N. gonorrhoeae* due to its high feasibility.

Although the RPA-Cas12a system was successfully applied to *N. gonorrhoeae* detection, our study has several limitations. According to the TwistDx instruction manual, primer amplification performance cannot be judged only by primer sequence, candidate primers need to be tested and screened. However, we screened only six candidate primer pairs and prioritized primer specificity over sensitivity during primer screening. The RPA-Cas12a-based *N. gonorrhoeae* detection system developed in this study was less sensitive than the previously reported 1 copy/ μ L, possibly due to the RPA primer selection. Encouragingly, there were no false

negatives when compared with the results of the gold standard traditional culture method, indicating that the sensitivity of 5 pg/ μ L of this method is adequate for clinical detection. Constant shaking during the RPA reaction has been shown to enhance the RPA amplification rate and improve sensitivity (Li et al. 2018), which could help to improve the LoD of the RPA-Cas12a system. To improve the sensitivity of the RPA-Cas12a system, further research should focus on primers with higher amplification efficiency and RPA reactions with constant shaking. In this study, only species identification of *N. gonorrhoeae* was carried out. With the high levels of antimicrobial resistance (AMR) (Unemo et al. 2019), further development of a test based on Cas12a that can simultaneously detect AMR, will make a great




Fig. 4 Clinical validation. **a** Clinical validation of the RPA-Cas12a-fluorescent assay. 24 clinical samples from individuals suspected of having gonorrhoea were examined by the RPA-Cas12a-fluorescent assay. AU, arbitrary unit. **b** Clinical validation of the RPA-Cas12a-LFD assay. 24 clinical samples from individuals suspected of having gonorrhoea were examined by the RPA-Cas12a-LFD assay

Table 3 Clinical validation of the RPA- Cas12a system

	Clinical samples																							
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
RPA-Cas12a fluorescent assay	+ ^a	+	+	- ^b	+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
RPA-Cas12a LFD assay	+	+	+	-	+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-
Traditional culture	+	+	+	-	+	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	-

^a *N. gonorrhoeae* was detected.

^b *N. gonorrhoeae* was not detected

contribution to both antibiotic selection and treatment of multidrug-resistant gonorrhoea.

In conclusion, we demonstrated a promising point-of-care test for *N. gonorrhoeae* with high sensitivity and specificity based on the RPA-Cas12a detection system. Moreover, the RPA-Cas12a-LFD assay can be widely available, which may contribute significantly to *N. gonorrhoeae* management in settings with limited medical facilities.

Supplementary Information

The online version contains supplementary material available at <https://doi.org/10.1186/s13568-023-01554-7>.

Additional file 1: Table S1. The NCBI accession numbers of *porA* gene of *N. meningitidis*. **Table S2.** Primers used in this study. **Table S3.** crRNA sequence (5-3'). **Table S4.** ssDNA reporter sequence.

Acknowledgements

Thanks to Professor Xiaobing Zhang for donating *N. gonorrhoeae* reference strain (ATCC 43069).

Author contributions

QT, PY and SH conceived and designed the research. QT, PY, LX, and CL conducted experiments and analyzed data. QT and XC contributed new reagents or analytical tools. QT and SH wrote the manuscript. All authors read and approved the manuscript.

Funding

This work was sponsored by the National Natural Science Foundation of China (Grant No. 82072346 and No. 82072349) and the Natural Science Foundation of Chongqing, China (No. cstc2020cyj-msxmX0519).

Availability of data and materials

Data available on request from the authors.

Declarations**Ethics approval and consent to participate**

The study was approved by the Institutional Medical Ethics Review Board of the First Affiliated Hospital of Chongqing Medical University. Isolates and clinical samples were collected at the microbiology laboratory of the hospital, with no contact with the patients. Therefore, the Ethics Review Board determined that informed consent was not required.

Consent for publication

All authors agree to be published.

Competing interests

The authors declared no competing interests.

Author details

¹Department of Clinical Laboratory Medicine, The First Affiliated Hospital of Chongqing Medical University, No. 1 Friendship Road, Yuzhong District, Chongqing 400016, People's Republic of China. ²Department of Burn and Plastic Surgery, The First Affiliated Hospital of Chongqing Medical University, No. 1 Friendship Road, Yuzhong District, Chongqing 400016, People's Republic of China. ³Department of Clinical Laboratory Medicine, Chongqing Shapingba District Chenjiaqiao Hospital, Chongqing 401331, People's Republic of China.

Received: 25 April 2023 Accepted: 1 May 2023

Published online: 27 May 2023

References

- Chen JS, Ma E, Harrington LB, Da Costa M, Tian X, Palefsky JM, Doudna JA (2018) CRISPR-Cas12a target binding unleashes indiscriminate single-stranded DNase activity. *Science* 360(6387):436–439. <https://doi.org/10.1126/science.aar6245>
- Gootenberg JS, Abudayyeh OO, Lee JW, Essletzbichler P, Dy AJ, Joung J, Verdine V, Donghia N, Daringer NM, Freije CA, Myhrvold C, Bhattacharyya RP, Livny J, Regev A, Koonin EV, Hung DT, Sabeti PC, Collins JJ, Zhang F (2017) Nucleic acid detection with CRISPR-Cas13a/C2c2. *Science* 356(6336):438–442. <https://doi.org/10.1126/science.aam9321>
- Harrington LB, Burstein D, Chen JS, Paez-Espino D, Ma E, Witte IP, Cofsky JC, Kyrpides NC, Banfield JF, Doudna JA (2018) Programmed DNA destruction by miniature CRISPR-Cas14 enzymes. *Science* 362(6416):839–842. <https://doi.org/10.1126/science.aav4294>
- Hjelmevoll SO, Olsen ME, Sollid JU, Haaheim H, Unemo M, Skogen V (2006) A fast real-time polymerase chain reaction method for sensitive and specific detection of the *Neisseria gonorrhoeae* *porA* pseudogene. *J Mol Diagn* 8(5):574–581. <https://doi.org/10.2353/jmoldx.2006.060024>
- Kersting S, Rausch V, Bier FF, von Nickisch-Roseneck M (2014) Multiplex isothermal solid-phase recombinase polymerase amplification for the specific and fast DNA-based detection of three bacterial pathogens. *Mikrochim Acta* 181(13–14):1715–1723. <https://doi.org/10.1007/s00604-014-1198-5>
- Li J, Macdonald J, von Stetten F (2018) Review: a comprehensive summary of a decade development of the recombinase polymerase amplification. *Analyst* 144(1):31–67. <https://doi.org/10.1039/c8an01621f>
- Lovett A, Duncan JA (2018) Human immune responses and the natural history of *Neisseria gonorrhoeae* infection. *Front Immunol* 9:3187. <https://doi.org/10.3389/fimmu.2018.03187>
- Mangold KA, Regner M, Tajuddin AM, Jennings L, Du H, Kaul KL (2007) *Neisseria* species identification assay for the confirmation of *Neisseria gonorrhoeae*-positive results of the COBAS Amplicor PCR. *J Clin Microbiol* 45(5):1403–1409. <https://doi.org/10.1128/jcm.00834-06>
- Qin X, Zhao Y, Chen W, Wu X, Tang S, Li G, Yuqi Y, Cao W, Liu X, Huang J, Yang J, Chen W, Tang W, Zheng H (2019) Changing antimicrobial susceptibility and molecular characterisation of *Neisseria gonorrhoeae* isolates in Guangdong, China: in a background of rapidly rising epidemic. *Int J Antimicrob Agents* 54(6):757–765. <https://doi.org/10.1016/j.ijantimicag.2019.08.015>
- Rowley J, Vander Hoorn S, Korenromp E, Low N, Unemo M, Abu-Raddad LJ, Chico RM, Smolak A, Newman L, Gottlieb S, Thwin SS, Broutet N (2019) *Chlamydia, gonorrhoea, trichomoniasis and syphilis*: global prevalence and incidence estimates, 2016. *Bull World Health Organ* 97(8):548–562p. <https://doi.org/10.2471/blt.18.228486>
- Stenger MR, Pathela P, Schumacher C, Burghardt N, Amiya R, Madera R, Nguyen TQ, Torrone E (2021) Trends in HIV prevalence by self-report among MSM diagnosed and reported with gonorrhoea in six United States jurisdictions from 2010 to 2019. *AIDS* 35(15):2523–2530. <https://doi.org/10.1097/qad.0000000000003067>
- Stevens JS, Criss AK (2018) Pathogenesis of *Neisseria gonorrhoeae* in the female reproductive tract: neutrophilic host response, sustained infection, and clinical sequelae. *Curr Opin Hematol* 25(1):13–21. <https://doi.org/10.1097/moh.0000000000000394>
- Su WH, Tsou TS, Chen CS, Ho TY, Lee WL, Yu YY, Chen TJ, Tan CH, Wang PH (2011) Are we satisfied with the tools for the diagnosis of gonococcal infection in females? *J Chin Med Assoc* 74(10):430–434. <https://doi.org/10.1016/j.jcma.2011.08.012>
- Swarts DC (2019) Making the cut(s): how Cas12a cleaves target and non-target DNA. *Biochem Soc Trans* 47(5):1499–1510. <https://doi.org/10.1042/bst20190564>
- Unemo M, Seifert HS, Hook EW, Hawkes S, Ndowa F, Dillon JR (2019) Gonorrhoea. *Nat Rev Dis Primers* 5:1–79. <https://doi.org/10.1038/s41572-019-0128-6>
- Vallye LM, Egli-Gany D, Wand H, Pomat WS, Homer CSE, Guy R, Silver B, Rumbold AR, Kaldor JM, Vallye AJ, Low N (2021) Adverse pregnancy and neonatal outcomes associated with *Neisseria gonorrhoeae*: systematic review and meta-analysis. *Sex Transm Infect* 97(2):104–111. <https://doi.org/10.1136/sextrans-2020-054653>
- Vigué L, Eyre-Walker A (2019) The comparative population genetics of *Neisseria meningitidis* and *Neisseria gonorrhoeae*. *PeerJ* 7:e7216. <https://doi.org/10.7717/peerj.7216>
- WHO (2021) Global progress report on HIV, viral hepatitis and sexually transmitted infections, 2021. World Health Organization. <https://www.who.int/publications/i/item/9789240027077>. Accessed 15 July 2021
- WHO (2016) WHO guidelines for the treatment of *Neisseria gonorrhoeae*. World Health Organization. <https://www.who.int/publications/i/item/9789241549691>. Accessed 1 January 2016
- Xiong D, Dai W, Gong J, Li G, Liu N, Wu W, Pan J, Chen C, Jiao Y, Deng H, Ye J, Zhang X, Huang H, Li Q, Xue L, Zhang X, Tang G (2020) Rapid detection of SARS-CoV-2 with CRISPR-Cas12a. *PLoS Biol* 18(12):e3000978. <https://doi.org/10.1371/journal.pbio.3000978>
- Yin YP, Han Y, Dai XQ, Zheng HP, Chen SC, Zhu BY, Yong G, Zhong N, Hu LH, Cao WL, Zheng ZJ, Wang F, Zhi Q, Zhu XY, Chen XS (2018) Susceptibility of *Neisseria gonorrhoeae* to azithromycin and ceftriaxone in China: a retrospective study of national surveillance data from 2013 to 2016. *PLoS Med* 15(2):e1002499. <https://doi.org/10.1371/journal.pmed.1002499>
- Zetsche B, Gootenberg JS, Abudayyeh OO, Slaymaker IM, Makarova KS, Essletzbichler P, Volz SE, Joung J, van der Oost J, Regev A, Koonin EV, Zhang F (2015) Cpf1 is a single RNA-guided endonuclease of a class 2 CRISPR-Cas system. *Cell* 163(3):759–771. <https://doi.org/10.1016/j.cell.2015.09.038>

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.