

POSTER PRESENTATION

Open Access

# Left atrial wall imaging using a novel black-blood late gadolinium enhancement sequence

Maxine Tang<sup>1,2\*</sup>, Tamer Basha<sup>1</sup>, Sophie Berg<sup>1</sup>, Kraig V Kissinger<sup>1</sup>, Beth Goddu<sup>1</sup>, Warren J Manning<sup>1,2</sup>, Reza Nezafat<sup>1</sup>


From 19th Annual SCMR Scientific Sessions  
Los Angeles, CA, USA. 27-30 January 2016

## Background

Late gadolinium enhancement (LGE) is commonly used clinically to identify myocardial infarct or scar as a hyper-enhanced region. However, the blood pool also remains bright in such images, masking adjacent bright tissue. LGE has also been used extensively to image fibrosis in the left atrium (LA) in patients with atrial fibrillation. However, identification of LA wall in LGE is challenging due to the thin wall of LA and high signal of the blood pool in the LA. We recently developed a novel black-blood LGE (BB-LGE) sequence (Fig. 1) which uses an optimized combination of inversion-recovery and  $T_2$ -preparation ( $T_2$ -prep) prepulses to simultaneously null the healthy myocardium and blood signals by leveraging the higher  $T_2$  of blood. In this study, we sought to assess the feasibility of BB-LGE in imaging of LA.

## Methods

Seven patients (6 male, 65+/-17 years) with a history of atrial fibrillation were recruited to participate in an IRB-approved study to compare BB-LGE vs. LGE for imaging of the atrial wall. All images were acquired at 1.5T (Philips Achieva). LGE and BB-LGE scans were acquired 15-20 minutes after injection of 0.1 mmol/kg of Multi-Hance contrast agent. Both sequences were acquired with identical imaging parameters and spatial resolution; BB-LGE required an additional  $T_2$ -prep to null the blood signal. The parameters of the BB-LGE sequence ( $\Delta t_1$ ,  $\Delta t_2$ ,  $\Delta t_3$  in Fig. 1) were determined in a short ~30 sec navigator scan by a technologist prior to the BB-LGE scan. Images were subjectively scored by an observer according to percentage of the LA border that was clearly defined. A score of 1 was given for scans in which greater than 75% of the border was well-defined.



A score of 2 was given for 50-75%, a score of 3 was given for 25-50%, and a score of 4 was given for less than 25% of the LA border. The scores for each sequence were compared using a paired T-test analysis.

## Results

Figure 2 shows examples of BB-LGE and LGE images of the LA. The LA wall could be easily delineated in BB-LGE images. There was significant difference in the scores of LA border definition for the BB-LGE images vs. the LGE images ( $p = 0.02$ ). The BB-LGE image scores were lower than the LGE image scores (1.38+/-0.74 vs. 2.25+/-0.89), indicating that a higher percentage of the LA border was well-defined on the BB-LGE images.

## Conclusions


Our proposed BB-LGE sequence has the potential to provide improved visualization of the LA wall compared to LGE protocols due to its ability to simultaneously null both the healthy myocardium and blood pool signals. This preliminary data demonstrates the need for further study of this sequence in patients with suspected LA fibrosis.

## Authors' details

<sup>1</sup>Cardiology, Beth Israel Deaconess Medical Center, Boston, MA, USA.

<sup>2</sup>Radiology, Beth Israel Deaconess Medical Center, Boston, MA, USA.

<sup>1</sup>Cardiology, Beth Israel Deaconess Medical Center, Boston, MA, USA  
Full list of author information is available at the end of the article



Published: 27 January 2016

doi:10.1186/1532-429X-18-S1-P62

Cite this article as: Tang et al.: Left atrial wall imaging using a novel black-blood late gadolinium enhancement sequence. *Journal of Cardiovascular Magnetic Resonance* 2016 **18**(Suppl 1):P62.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at  
[www.biomedcentral.com/submit](http://www.biomedcentral.com/submit)

 BioMed Central