

POSTER PRESENTATION

Open Access

# A compressed sensing multi-slice cine CMR approach for the accurate assessment of left ventricular volumes and function

Gabriella Vincenti<sup>1\*</sup>, Davide Piccini<sup>2,3</sup>, Pierre Monney<sup>1</sup>, Jerome Chaptinel<sup>3</sup>, Tobias Rutz<sup>1</sup>, Simone Coppo<sup>3</sup>, Michael O Zenge<sup>4</sup>, Michaela Schmidt<sup>4</sup>, Mariappan S Nadar<sup>5</sup>, Pascal Chevre<sup>1,6</sup>, Matthias Stuber<sup>3</sup>, Juerg Schwitler<sup>1</sup>

From 17th Annual SCMR Scientific Sessions  
New Orleans, LA, USA. 16-19 January 2014

## Background

CMR is generally accepted as the gold standard for left ventricular (LV) volumes and function assessment. The conventional CMR approach involves several breath-holds to cover the entire heart with short-axis acquisitions. Recently, compressed sensing (CS) techniques emerged as a means to considerably accelerate data acquisition. CS principally relies on: 1) transform sparsity, 2) incoherence of undersampling artifacts, and 3) nonlinear reconstruction. **PURPOSE:** To compare a novel CS-based single breath-hold multi-slice cine technique with the standard multi-breath-hold technique for the assessment of LV volumes and function.

## Methods

The novel single breathhold multi-slice CS cine sequence was implemented on a 1.5T MAGNETOM Aera (Siemens Healthcare, Germany). It realizes incoherent sampling by distributing the readouts pseudo-randomly on the Cartesian k space grid and by applying a pseudo-random offset from frame-to-frame during prospective ECG-triggering. Thus, 3 long-axis and 4 short-axis slices (Figure 1A) were acquired in a single breath-hold of 14 heart beats (temporal/spatial resolution: 30 ms/1.5 × 1.5 mm<sup>2</sup>, respectively, acceleration factor: 11.0). The CS cine data were analyzed by the Siemens Argus 4DVF software (Figure 1B) which is based on a 3D LV-model that takes the cyclic motion of the mitral valve plane (yellow plane in 2B) into account. For comparison, a conventional stack of cine SSFP images was acquired (temporal/spatial resolution 40 ms/1.2 × 1.6 mm<sup>2</sup>, slice thickness:

8 mm, gap: 2 mm) and analyzed by the Siemens Argus VF software. As a reference for the LV stroke volume (LVS<sub>V</sub>), the aortic flow (AoFlow) was measured by a phase-contrast acquisition (temporal/spatial resolution 40 ms/1.8 × 1.8 mm<sup>2</sup>).

## Results

Twelve volunteers (75% male, age 33 ± 7 y) and 14 patients (85% male, age 67 ± 10 y) with different LV pathologies were studied. Good image quality of the single breathhold multi-slice CS acquisitions was obtained in 10/12 volunteers and 12/14 patients and all CS data were of adequate quality for 4D-analysis. LVEF<sub>CS</sub> and LVEF<sub>Standard</sub> were similar (49.3 ± 16.4% vs 50.6 ± 16.4%, respectively, p = 0.17) with an excellent correlation (Figure 2). Correlation of LVS<sub>VCS</sub> with AoFlow was also excellent (r = 0.93, p < 0.0001). For LVS<sub>VStandard</sub> the correlation was r = 0.70 (p < 0.01). Bland-Altman analyses showed underestimation of CS vs standard by 10.5 ml and 7.2 ml for LV end-diastolic (p < 0.005) and LV stroke volume (p < 0.002), respectively, no difference for LV end-systolic volume (p = ns).

## Conclusions

This is the first clinical experience with a novel CS-based cine sequence. These preliminary results demonstrate the potential of this fast single breath-hold technique to replace the multi-breath-hold standard technique for the assessment of LV volumes and function.

## Funding

D. Piccini, M.Zenge, M.Schmidt and M. Nadar are Siemens employees. Other funding than that was not received.

<sup>1</sup>Division of Cardiology and Cardiac MR Center, University Hospital of Lausanne (CHUV), Lausanne, Switzerland  
Full list of author information is available at the end of the article




Figure 1




Figure 2

#### Authors' details

<sup>1</sup>Division of Cardiology and Cardiac MR Center, University Hospital of Lausanne (CHUV), Lausanne, Switzerland. <sup>2</sup>Advanced Clinical Imaging Technology, Siemens Healthcare IM BM PI, Lausanne, Switzerland. <sup>3</sup>Department of Radiology, University Hospital (CHUV) and University of Lausanne (UNIL)/Center for Biomedical Imaging (CIBM), Lausanne,

Switzerland. <sup>4</sup>MR Applications and Workflow Development, Healthcare Sector, Siemens AG, Erlangen, Germany. <sup>5</sup>Siemens Corporate Technology, Princeton, New Jersey, USA. <sup>6</sup>Department of Radiology, University Hospital Lausanne (CHUV), Lausanne, Switzerland.

Published: 16 January 2014

doi:10.1186/1532-429X-16-S1-P376

**Cite this article as:** Vincenti et al.: A compressed sensing multi-slice cine CMR approach for the accurate assessment of left ventricular volumes and function. *Journal of Cardiovascular Magnetic Resonance* 2014 16 (Suppl 1):P376.

**Submit your next manuscript to BioMed Central and take full advantage of:**

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at  
[www.biomedcentral.com/submit](http://www.biomedcentral.com/submit)

